
50 50

## Puntos clave. *Punti chiavi.*

## Points clé. *Key points.*

1

Conector: MP-5.  
*Connettore: MP-5.*  
Connecteur: MP-5.  
*Connector: MP-5.*

2

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# pH electrode 50 50

## Introduction

The **50 50** is a universal pH electrode to be used with the HACH portable pH-meter. These are low maintenance electrodes, as they do not require monitoring of electrolyte levels.

## Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it.

 **MP-5 connector (IP 67).** Connect and disconnect by pushing. Do not unscrew it.

## Calibration

All pH electrodes must be calibrated with buffer solutions before measuring.  See pH-meter instruction manual.

## Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases. In field measurements, these types of electrodes allow stirring the sample directly holding the electrode by the handle.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

## WARNING

- Do not leave the electrode immersed in distilled water.
- Do not bend the cable on the handle level.
- Never unscrew the MP-5 connector.

## Maintenance

- Maintain the membrane's surface ⑧ clean.
- When an electrode gets contaminated, the simply rinsing with distilled water is not sufficient. The most suitable cleaning solution is which acts in a most selective way over the contamination (see HACH regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in HACH cleaning electrode solution, code C20C370, or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrodes regenerated solutions.

## Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode, cable or connector.

**Substitute the electrode.**

... fortuitous reading.

Broken cable.

**Substitute the electrode.**

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm and/or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

**Clean, regenerate or replace the electrode.**

## Storage

Place the electrode in the storage protector ④ with KCl 3M sat. AgCl.

## Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not response appropriately after regeneration, must be replaced by a new one. To avoid accidental breakage of the glass body electrodes, the use of the **working protector** is recommended.

 *If an electrode has too short lifespan, less than three months, ask us about it. Probably there are more suitable electrode for the application.*

## Specifications

|                      |                                                                                           |
|----------------------|-------------------------------------------------------------------------------------------|
| pH measuring range   | 0 ... 14                                                                                  |
| Operating temp.      | 0 ... 80 °C                                                                               |
| Reference element    | Ag wire coated with AgCl                                                                  |
| Diaphragms           | ceramic                                                                                   |
| Electrolyte          | gel                                                                                       |
| Body material        | polycarbonate                                                                             |
| Min. immersion depth | 18 mm  |

## Applications:

- Aqueous samples in general.

## Limits

Low-conductivity, dirty or viscous solutions.

Products with colloids or solids in suspension.

Samples containing sulphides, reducing sugars or other substances which react with the silver ions.

Temperatures > 80 °C.

## Warranty

HACH guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

**Limitations:**

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

## CERTIFICATE OF QUALITY

The attached electrode **50 50** has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s

**HACH COMPANY World Headquarters**

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.

Tel. (970) 669-3050

(800) 227-4224 (U.S.A. only)

Fax (970) 669-2932

orders@hach.com

www.hach.com

**HACH LANGE GMBH**

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D-40549 Düsseldorf, Germany

Tel. +49 (0) 2 11 52 88-320

Fax +49 (0) 2 11 52 88-210

info@de.hach.com

www.de.hach.com

**HACH LANGE Srl**

6, route de Compois

1222 Vérenaz

SWITZERLAND

Tel. +41 22 594 6400

Fax +41 22 594 6499