



**POL-EKO**  
Perfect Environment

(+48) 32 453 91 70  
info@pol-eko.com.pl  
www.pol-eko.com.pl

**POL-EKO A. Polok-Kowalska sp.k.**  
44-300 Wodzisław Śląski  
ul. Kokoszycka 172 C

## Data sheet

## Emergency power supply ZA II 8H



The photo above is for reference only, may show additional options not included in standard equipment. The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.

### Standard equipment

- converter with the function of charging the batteries
- battery/ies
- castors with brake
- visual and acoustic signal informing about the operating status
- electric socket type E (230V)
- user manual in English



**POL-EKO**  
Perfect Environment

(+48) 32 453 91 70  
info@pol-eko.com.pl  
www.pol-eko.com.pl


**POL-EKO A. Polok-Kowalska sp.k.**  
44-300 Wodzisław Śląski  
ul. Kokoszycka 172 C

## TECHNICAL DATA

display	LED
housing	powder coated sheet

### overall dims [mm]

width A	675
height B	670
depth C	665
weight [kg]	173



## BASIC PARAMETERS

operating temperature [°C]	+10...+30
is compatible with the model	CHL/ST 500-1450, ILW 240-750
back-up time* [h]	8
number of batteries [pcs.]	2

## ELECTRICAL PARAMETERS

voltage**	230 V 50 Hz
nominal power [W]	3 000
warranty	24 months
manufacturer	POL-EKO

\* the approximate time of maintaining the operation of the unit with the emergency power supply depends on the parameters of the environment, the degree of filling the chamber, etc.